



Needle-Free IVF: FAQ

DOES IT HURT TO APPLY THE NEEDLE-FREE DEVICE?

INSERTION OF THE DEVICE IS QUICK AND SIGNIFICANTLY LESS PAINFUL THAN A TYPICAL SUBCUTANEOUS INJECTION. THIS DEVICE REDUCES THE NUMBER OF NEEDLE INJECTIONS YOU WILL NEED DURING AN IVF CYCLE

WHAT MEDICATIONS CAN BE USED WITH NEEDLE-FREE IVF?

ALL MEDICATIONS THAT ARE PRESCRIBED TO BE GIVEN SUBCUTANEOUSLY. THIS DEVICE REQUIRES A LUER LOCK (TWIST ON) SYRINGE. FOR MEDICATIONS IN A PEN OR WITH A NON-REMOVABLE NEEDLE, AN ADDITIONAL STEP MAY BE REQUIRED.

INTRAMUSCULAR MEDICATIONS ARE NOT COMPATIBLE WITH NEEDLE-FREE IVF

HOW LONG DOES THE DEVICE LAST?

THE DEVICE LASTS FOR 9-12 DAYS AND CAN BE USED FOR THE DURATION OF YOUR STIMULATION. TYPICALLY, THE DEVICE WILL BE REPLACED ONCE DURING YOUR STIM.

WHAT ARE MY RESTRICTIONS?

YOU MAY CONTINUE ALL OF YOUR NORMAL LIFESTYLE ACTIVITIES. HOWEVER, WE RECOMMEND USING CAUTION WITH SWIMMING OR HIGH-SWEAT ACTIVITIES SUCH AS SAUNAS, STEAM ROOMS, OR HOT YOGA, AS THESE CAN WEAKEN THE DRESSING

HOW DO I MAKE SURE I ADMINISTERED MY MEDICATION CORRECTLY?

OUR STAFF IS AVAILABLE TO APPLY THE DEVICE FOR YOUR FIRST USE. WE WILL GUIDE YOU THROUGH MEDICATION ADMINISTRATION WITH VIDEOS AND TOOLS TO MAKE THE PROCESS SEAMLESS. A SALINE FLUSH WILL BE PROVIDED TO FLUSH THE TUBING AFTER EACH MEDICATION ADMINISTRATION